

BLUSTER 3 (ST6 3 wheel)

Item No. 950100600

drive

DeVilbiss
HEALTHCARE



- Luxury 3 wheel **15 km/h** scooter
- Pneumatic tires and suspension all around
- Fully adjustable Captain chair
- Available with seat suspension
- Manoeuvrable and comfortable
- Infinitely adjustable steering column
- Delta steering

Technical Data

Order number	950100600
Colour	anthracite
Overall width	72 cm
Overall length	158 cm
Overall height	136 cm
Seat width	58 cm
Seat depth	56 cm
Seat height	48 / 50.5 - 53 cm
Width between armrest	55 cm
Speed	15 km/h
Battery capacity 2x	80 Ah

Motor capacity	800 W
Continuous-power	800 W
Peak power	1600 W
Gradeability	16 %
Clearance	13 cm
Front tires	3.5 - 8
Rear tires	3.5 - 8
Battery Voltage	80 V
Weight without battery	107 kg
Max. Load	180 kg