

# Hybrid power complete system with pump

Item No. HYB/POWER



The Hybrid-Power mattress combines the pressure redistribution properties of high specification foam, with alternating pressure relief.

- Cover can withstand a higher concentration of cleaning substances and offers greater resistance to abrasion
- High frequency welded seams protect against fluid ingress
- 360° zip to prolong longevity • High specification foam combined with alternating cells
- Zonal mattress surface gives additional support to vulnerable
- NEW full width heel zone areas
- Recommend as therapeutic support for all decubitus categories (1 to 4\*) and for prophylaxis in patients up to a very high decubitus risk

Important

\*For people who can change their lying position regularly themselves or with help!

## Technical Data

Order number	HYB/POWER
Overall width	88 cm
Overall length	198 cm

Overall height	16 cm
User weight	254 kg
Max. Load	254 kg